## **Guide to Selecting the Appropriate Rheometer and Spindle**

In order to determine which rheometer is best for the sample, the user must first know the viscosity range at which the sample falls under. If the viscosity falls in the range between 15-2,000,000 cPs, the user should select the DV3TLV rheometer. Some common fluids that fall in this viscosity range include inks, oils and solvents. If the viscosity falls in the range between 2,000,000-80,000,00 cPs, the user should select the DV3THA rheometer. Some common fluids that fall in this viscosity range include gels, chocolates, and epoxies. After selecting the appropriate rheometer, the user should then use the table below to select the appropriate spindle.

	DV3TLV				DV3THA		
	Spindles	Viscosity Range [cP]	Sample Size [mL]		Spindles	Viscosity Range [cP]	Sample Size [mL]
Small Sample Adapter	SC4-21	2.4-46,865	7.1	Small Sample Adapter	SC4-21	2.4-46,865	7.1
	SC4-27	46,865-234,325	10.4		SC4-27	46,865- 234,325	10.4
ΓΛ	LV1	15-20,000		11م	LV1	15-20,000	
	LV2	20,000-100,000			LV2	20,000- 100,000	
	LV3	100,000-400,000			LV3	100,000- 400,000	
	LV4	400,000- 2,000,000			LV4	400,000- 2,000,000	
НА	HA1	200-40,000		НА	HA1	200-40,000	
	HA2	40,000- 160,000			HA2	40,000- 160,000	
	HA3	160,000-400,000			HA3	160,000- 400,000	
	HA4	400,000-800,000			HA4	400,000- 800,000	
	HA5	800,000- 1,600,000			HA5	800,000- 1,600,000	
	HA6	1,600,000- 2,000,000			HA6	1,600,000- 4,000,000	
T-Bar (Helipath)	TA	156-62,500		T-Bar (Helipath)	TA	4,000-8,000	
	ТВ	62,500-124,800			ТВ	8,000-20,000	
	TC	124,800-312,000			TC	20,000-40,000	
	TD	312,000-624,000			TD	40,000- 8,000,000	
	TE	624,000- 1,500,000			TE	8,000,000- 20,000,000	
	TF	1,500,000- 2,000,000			TF	20,000,000- 40,000,000	